



Propel Half Label Comparison - Vial Label

PROPEL RX



PHARMACLIK RX



- 1 Prescriber the prescriber's salutation and name have been moved under the prescription number (authorization number in British Columbia).
- *NEW* DIS ID In Newfoundland and Labrador, this is the Create Rx ID from the Rx Detail EHR tab. If the network is down when the label is printed, this field is blank. This does not print for shrink labels.
- 3 Fill Date the Fill Date has been moved to below the patient's name.
- **Drug Name** Trade and Trade & Innovator options for the drug name have been removed to adhere to labelling standards. The drug name continues to print according to the settings in the Drug Folder, Group Folder, and Rx Detail. Mixture prescriptions will always display the generic and trade name. The generic name is replaced with TallMan if it exists in the Drug Folder (not applicable for mixtures).
 - In New Brunswick, the MFR also prints beside the strength. In Nova Scotia, for insulin prescriptions, the trade name always prints on the first line and the generic name on the second.
- Quantity the Quantity has been moved to the same row as the dosage form, manufacturer, and DIN. For Owe Balance labels, this is the amount that's being dispensed for the Owe Balance. For methadone prescriptions, this is the number of millilitres per dose up to 1 decimal place.





Last Updated: September 2023

- Refills Expire/No Refills After in Nova Scotia and Newfoundland and Labrador, the Rx expiry date must print on the label; the Rx Life Expiry preference is defaulted to Written Rx Date and disabled, and the Print Rx Expiry preference must be set to either # Refill(s) Expire MMM DD, YYYY or No Refills After MMM DD, YYYY. In all other provinces, the Print Rx Expiry preference can be set to None in addition to the aforementioned options.
- **Qty Left/Rep** Qty Left prints instead of Rep if the **Use Qty Left for Repeats** preference is ON. Qty Left always prints if a partial fill is left. In Manitoba, Qty Left always prints.
- 8 SIG the SIG has been moved to the bottom of the vial label. The maximum number of characters remains the same (193 mixed or 187 uppercase). The text shrinks as needed to fit within 3 lines.
- *NEW* Barcode a prescription barcode prints vertically on the right side of the label. Accucount barcodes now print on an auxiliary label instead of the vial label.
- *NEW* Label Indicator the label number and the total number of labels are indicated in the bottom right corner if a split is specified in Rx Detail or Packaging (format: label #/total # of labels). For owe balances, the label indicator reflects the number of owe balance labels printed. If there is no split, this field does not appear. This does not print for methadone, Additional Instructions, supplementary, immunization, or Hold labels. This field can be updated with an Amend or Amend Next.





Propel Daily Laser Comparison - Vial Label

PROPEL RX



PHARMACLIK RX



- 1 Drug Name Trade and Trade & Innovator options for the drug name have been removed to adhere to labelling standards. The drug name continues to print according to the settings in the Drug Folder, Group Folder, and Rx Detail. Mixture prescriptions will always display the generic and trade name. The generic name is replaced with TallMan if it exists in the Drug Folder (not applicable for mixtures).
 - In New Brunswick, the MFR also prints beside the strength. In Nova Scotia, for insulin prescriptions, the trade name always prints on the first line and the generic name on the second.
- *NEW* Area indicates the name of the area that the patient belongs to in the group.
- **Qty Left/Rep** Qty Left prints instead of Rep if the **Use Qty Left for Repeats** preference is ON in Half Label Preferences. Qty Left always prints if a partial fill is left. In Manitoba, Qty Left always prints.
- *NEW* Refills Expire/No Refills After in Nova Scotia and Newfoundland and Labrador, the Rx expiry date must print on the label; the Rx Life Expiry preference is defaulted to Written Rx Date and disabled, and the Print Rx Expiry preference must be set to either # Refill(s)

 Expire MMM DD, YYYY or No Refills After MMM DD, YYYY. In all other provinces, the Rx expiry date does not print on the Daily Laser vial label.
- *NEW* Total Paid in Ontario, the price of the prescription appears here.

